

PROJECT AT KUMC PRECISION NEURAL  
DYNAMICS LAB:  
**OBJECT TRACK MATRIX**

---

Presenter: Thresa Kelly

Authors: Thresa Kelly, Ryan Coppens, Dr. Alice Bean, & Dr. Adam Rouse

# INTRODUCTION

---

# Project Goal

- Design an activity board with software for monkey reach-and-grasp experiments
  - 4x3 grid of motion sensors (12 total)
  - LED lights up object → indicate to grab object
  - Software to customize experiment and record real-time data

# Background

- People with neurological diseases, movement disorders, or amputations may have difficulty controlling their body
  - New technologies and robotics can help!
- University of Kansas (KU) **Department of Neurosurgery Precision Neural Dynamics Lab** → research on how the brain controls the body
  - Study brain/spinal cord activity in animals → better understand how groups of neurons coordinate for...
    - Arm reaching
    - Hand grasping
  - Develop computational methods for analyzing large datasets from experiments



# Project Motivation

- Object track matrix → reach & grasp experiments with rhesus monkeys
- **Research:**
  - Study neural encoding of precision movements → target-independent and target-dependent neural signal limitations of speed and precision
    - Independent = neurons that fire for any movement in a given motor area
    - Dependent = level of activity of specific neurons for movement direction
  - Neural mapping of hand function → brain maps specific hand grasp, then sends signal to hand via spinal cord



# OVERALL DESIGN

---

# Requirements (1/2)

- **Physical board:**

- Objects must be able to be disconnected from the activity board.
- Rods must be able to support objects of different size, shape, and weight.
- The activity board must support 12 objects in a 4 x 3 grid.

- **Motion sensing:**

- Accelerometer must be able to detect motion in one dimension.
- Noise from the accelerometer must not be interpreted as an interaction from the monkey.
- The motion sensor package must be wirelessly attached to the end of the rod with the object.

# Requirements (2/2)


- **Lighting:**

- Only one LED should turn on at a time.
- The LED should turn on for a duration set by the user.

- **Software:**

- The user must be able to plan an experiment by determining the LED sequence and timing durations.
- The software must read from the hardware in real time.
- The timing resolution of collected data must be accurate to 10 ms.
- The software must determine if a trial was successful (the correct object was manipulated) or a failure (no object or the incorrect object was manipulated).






# Control Flow

- Main components:
  - Software
  - Activity board
  - Motion Sensors
- Motion sensors communicate with the activity board via Wi-Fi
- The activity board communicates with the software via serial communication via USB port


# MY CONTRIBUTIONS

---

# User Interface: Welcome



# User Interface: Experiment Details



The screenshot shows a software window titled "Object Track Matrix" with standard window controls (minimize, maximize, close). The window is divided into two main sections: "Experiment Details" and "Experiment Setup".

**Experiment Details**


- User Information**
  - Experiment title:
  - Experimenter name:
  - Subject name:
  - Date:  (with a date picker icon)
  - Notes:
  - Directory path:

**Experiment Setup**

- Port:  (dropdown menu)
- Time between trials:  (spinner) seconds
- Trial timeout:  (spinner) seconds
- Sample rate:  (spinner) Hz
- Sequence:
- Number of trials:  (spinner)

At the bottom center of the window is a  button.

# User Interface: Experiment




- Experiment data is saved to three files:
- <Directory Path>/<Experiment Name>/
  - Experiment\_Information.csv
  - Experiment\_Data\_Stream.csv
  - Experiment\_Data\_Stream.csv

# Software Internal Design

Class Name	Description
<b>MainWindow</b>	Top-level class and UI to request and submit experiment information, perform an experiment, and save data to files.
<b>GetUserInfo</b>	UI with user input fields for the experiment title, experimenter name, subject name, date, notes, and directory path. This information is accessed using getter functions.
<b>ExperimentSetup</b>	UI with user input fields for the serial port, time between trials, trial timeout, sample rate, and trial sequence. This information is accessed using getter functions.
<b>Commands</b>	Handles building and interpreting serial packets. Also has functions to convert integers to hexadecimal strings, check the command parameter requirements, and get the number of bytes of the packet components.
<b>SerialControl</b>	Handles reading and writing command packets through a serial port.
<b>ExperimentFileControl</b>	Writes files to the experiment directory to store the experiment information, streaming data, and trial status as *.csv files.

# Class Structure



# Serial Communication Protocol

Bytes	0	1	2	3	4	5	6	7
Component	STX	Command	ID	Data				ETX

- 1 byte = 8 bits = 1 ASCII character
- Building a packet:
  - Decimal integer → hexadecimal string ( 15 → "F")
  - Hex character → ASCII code ("F" → 70)

Decimal	Hex	Char	Decimal	Hex	Char	Decimal	Hex	Char	Decimal	Hex	Char
0	0	[NULL]	32	20	[SPACE]	64	40	@	96	60	`
1	1	[START OF HEADING]	33	21	!	65	41	A	97	61	a
2	2	[START OF TEXT]	34	22	"	66	42	B	98	62	b
3	3	[END OF TEXT]	35	23	#	67	43	C	99	63	c
4	4	[END OF TRANSMISSION]	36	24	\$	68	44	D	100	64	d
5	5	[ENQUIRY]	37	25	%	69	45	E	101	65	e
6	6	[ACKNOWLEDGE]	38	26	&	70	46	F	102	66	f
7	7	[BELL]	39	27	'	71	47	G	103	67	g
8	8	[BACKSPACE]	40	28	(	72	48	H	104	68	h
9	9	[HORIZONTAL TAB]	41	29	)	73	49	I	105	69	i
10	A	[LINE FEED]	42	2A	*	74	4A	J	106	6A	j
11	B	[VERTICAL TAB]	43	2B	+	75	4B	K	107	6B	k
12	C	[FORM FEED]	44	2C	,	76	4C	L	108	6C	l
13	D	[CARRIAGE RETURN]	45	2D	-	77	4D	M	109	6D	m
14	E	[SHIFT OUT]	46	2E	.	78	4E	N	110	6E	n
15	F	[SHIFT IN]	47	2F	/	79	4F	O	111	6F	o
16	10	[DATA LINK ESCAPE]	48	30	0	80	50	P	112	70	p
17	11	[DEVICE CONTROL 1]	49	31	1	81	51	Q	113	71	q
18	12	[DEVICE CONTROL 2]	50	32	2	82	52	R	114	72	r
19	13	[DEVICE CONTROL 3]	51	33	3	83	53	S	115	73	s
20	14	[DEVICE CONTROL 4]	52	34	4	84	54	T	116	74	t
21	15	[NEGATIVE ACKNOWLEDGE]	53	35	5	85	55	U	117	75	u
22	16	[SYNCHRONOUS IDLE]	54	36	6	86	56	V	118	76	v
23	17	[END OF TRANS. BLOCK]	55	37	7	87	57	W	119	77	w
24	18	[CANCEL]	56	38	8	88	58	X	120	78	x
25	19	[END OF MEDIUM]	57	39	9	89	59	Y	121	79	y
26	1A	[SUBSTITUTE]	58	3A	:	90	5A	Z	122	7A	z
27	1B	[ESCAPE]	59	3B	;	91	5B	[	123	7B	{
28	1C	[FILE SEPARATOR]	60	3C	<	92	5C	\	124	7C	
29	1D	[GROUP SEPARATOR]	61	3D	=	93	5D	]	125	7D	}
30	1E	[RECORD SEPARATOR]	62	3E	>	94	5E	^	126	7E	~
31	1F	[UNIT SEPARATOR]	63	3F	?	95	5F	_	127	7F	[DEL]




# Commands

Name	Number	Software		Hardware		Description
Testing		ID	Data			
PING	0					call and response to test communication
TEST LED	1	Object ID	Time (ms)			Command that turns on all LEDs to test functionality
BATTERY	2	Object ID		Object ID	Battery percent	Get the percent of the battery life of one motion sensor
Setup						
CALIBRATE	3	Object ID				Sends the object ID to be calibrated. The respective motion sensor is moved, and its ID saved.
NTRIALS	4		Number of trials			number of trials in the experiment
TRIAL	5	Object ID	Trial number			Sets the object ID for the trial. several of these commands comprise of the experiment object sequence.
SEPARATION	6		Time (ms)			Sets the time between trials. This is the time between an LED off and on.
TIMEOUT	7		Time (ms)			Sets the maximum time an LED can be on for a trial
SAMPLE RATE	8		Frequency (Hz)			Sets the sample rate
Experiment						
STREAM	9		True / False	Object ID	Trial number	Start (true) or stop (false) the experiment. Sends the object ID (1-12) of the object that was moved and the trial number. Object ID is 0 if no object is moved

# TESTING

---



# Read/Write Delay

- Computer time for:
  - the software to build a packet and send it via USB,  $t$
  - then for the microcontroller to interpret, execute, rebuild a packet, and write back to software,
  - Then for software to unpack and check the packet
- Goal:  $t < 10\text{ms}$
- **Average:  $4 \pm 1\text{ ms}$**

Trial	Did the correct LED ring turn on?	Did the motion sensor activate when moved?
1	Yes	Yes
2	Yes	Yes
3	Yes	Yes
4	Yes	Yes
5	Yes	Yes
6	Yes	Yes
7	Yes	Yes
8	Yes	Yes
9	Yes	Yes
10	Yes	Yes
11	Yes	Yes
12	No	No

# Test All Objects

- Sequence: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12.
- Write sequence to the activity board, then run experiment.
- Move each object when the LEDs light up.
- Goal: all objects must light up and respond correctly when moved.
- Result: all objects by #12 work.
  - Code bug when converting 12 to a hexadecimal character "C"...

# CONCLUSION

---

# Future work

- Fix bug with #12 integer to hex character conversion.
- Implement unused commands: TEST LED, BATTERY, and CALIBRATE
- Refactor serial communication protocol to include packet number
  - Protect against data loss
- Add real-time data visualization while running experiment

# QUESTIONS?

---

Thanks!

# DEMONSTRATION

---